

Carolina Chemistries Introduces Veintector

Carolina Liquid Chemistries introduces the Veintector, which is appropriate for doctor's offices or laboratories. This device can solve the technician's problem when searching for an "invisible" vein. It has an LED-absorbing hemoglobin property to locate and detect veins more easily than with the naked eye. The Veintector's slender shape, along with the tie belt, makes it easy to handle and battery-powered operation allows for easy carrying and storage.



Carolina Liquid Chemistries

1-877-722-8910

contactsales@carolinachemistries.com

www.carolinachemistries.com