



AIA[®]-1800

AUTOMATED IMMUNOASSAY ANALYZER

Fast, Flexible and Feature Rich

- ❖ Throughput of 180 results per hour
- ❖ Random access/continuous loading
- ❖ User defined programming
- ❖ Automated daily maintenance and record keeping
- ❖ Automatic inventory management of reagents and supplies
- ❖ True auto dilution from 2 to 625
- ❖ Sample clot detection
- ❖ Automated pretreatment
- ❖ Utilizes Tosoh's patented Unit Dose Test Cup technology - the same for all Tosoh AIA Systems



10-MINUTE (ST) and 40-MINUTE AIA-PACKS

Thyroid

- ST TSH
- ST TT3
- ST T4
- ST FT3
- ST FT4
- ST T-Uptake
- TSH 3rd Gen
- TPO Ab
- Tg Ab

Cardiac Markers

- ST CKMB
- ST Troponin I 2nd Gen
- ST Myoglobin

Reproductive

- ST β HCG
- ST Estradiol
- ST FSH
- ST LH II
- ST Prolactin
- ST Progesterone
- ST Testosterone

Tumor Markers

- ST AFP
- ST CEA
- ST PSA
- ST CA 125
- ST 27.29
- ST CA 19-9

Anemia

- ST Ferritin
- Vitamin B12
- Folate

Metabolic

- ST Cortisol
- ST C-Peptide
- ST GHG
- ST Insulin
- ST iPTH

Other

- ST β -2 Microglobulin
- ST IgE II
- ST PAP

TEST MENU

