

VITAMIN D

For Chemistry Analyzers

FDA-Cleared.

Recognizes both 25-OH Vitamin D2 and D3 equally.

This is a 2-part (R1 and R2) liquid ready-to-use reagent.

No sample pre-treatment required.

Good correlation with LC-MS/MS.

Traceability to NIST (SRM 972) and recovers DEQAS.

BENEFITS

- Easy: Runs on your chemistry analyzer
- Fast: Up to 800 tests/hour
- Economical: Far less expensive than IA analyzers

Description: R1 - 40mL / R2 - 10mL Cat.#: CCDZ888A-K

CALL TOLL FREE 877.722.8910

carolinachemistries.com



THE REAGENT COMPANY

Carolina Liquid Chemistries, Corp.
313 Gallimore Dairy Rd
Greensboro, NC 27409